

IMAGEWORKS EVA SELECT DIGITAL SENSOR SOLUTION



The *ImageWorks EVA Select Digital Sensor System*

has been completely redesigned for quality and performance

- Market-leading resolution (20+ line pairs per mm)
- Rounded corners for patient comfort
- Razor thin: 5.4 mm thickness
- Integrates with major software platforms
- USB interface works with your existing system




Gorgeous Images Plus a Great Price Equals Return on Investment



Universal Direct USB connection easily interfaces with your equipment

Need Imaging Software? We Have You Covered



EvaSoft Imaging software integrates all of your images into one system.

Easy and intuitive to use. No training needed. No mistakes.

Has capability to work with all major practice management platforms.

Powerful diagnostic capabilities.

One-click insurance claims allows you to work faster.

IMAGEWORKS EVA SELECT DIGITAL SENSOR SOLUTION

Technical Specifications:

	Size 2	Size 1
Gray shades:	4096 (12 bits)	4096 (12 bits)
Imager Resolution:	20+ lp/mm	20+ lp/mm
Optical Resolution:	26 lp/mm	26 lp/mm
Pixel Size:	19 microns x 19 microns	19 microns x 19 microns
Active Area:	36 mm (L) x 25.99 mm (W)	30.0 mm (L) x 20.0 mm (W)
Capsule Dimensions:	42.8 mm (L) x 30.5 mm (W) x 5.4 mm (Th)	36.7 mm (L) x 24.4 mm (W) x 5.4 mm (Th)
Image Size:	4.94 MB	3.16 MB



ImageWorks

Who we are:

We know that most dentists are also business owners who focus on getting the best value for their investment. Based in Elmsford, NY for over 35 years, ImageWorks has been manufacturing dental imaging solutions that provide great return on investment. This value stems not only from the quality of our products, but also the quality of our people.

The ImageWorks Full Dental Imaging End to End Solution



Imaging Software



Digital Intraoral Sensors



Intraoral X-Rays



Panoramic and
Cephalometric X-Rays



CBCT Systems