

2021 Products Catalog US & Canada

UPDATED 4/1/21





TABLE OF CONTENTS

ACTEON IMAGING

Intraoral Radiology	4
Extraoral Radiology	14
Intraoral Cameras	22

ACTEON EQUIPMENT

Ultrasonics	26
Air Polisher	43
Bliss Hand Instruments	48
Surgery	52
Curing Lights	
Curing Lights	სგ

ACTEON PHARMA

Gingival Retraction	71
Hemostasis & Tissue Regeneration	72
Hygiene	

ACTEON INTEGRATED

Ultrasonics	74
	79
	81

INTRAORAL RADIOLOGY 2D Imaging



X-MIND™ UNITY

X-MIND™ Unity is pre-wired for the integration of the SOPIX² digital sensor. Delivered with a long cone in standard version.

X-MIND™ UNITY - MOBILE STAND W0800140

X-MIND Unity Mobile Stand (XMU NOT INCLUDED)

X-MIND™ UNITY - TOP MOUNTING

XMU TOP 40 SHORT CONE

HXMU-TS

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 16"

XMU TOP 80 SHORT CONE

HXMU-TM

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 31"

XMU TOP 110 SHORT CONE

HXMU-TL

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Top Mount Extender Arm - 43"



INTRAORAL RADIOLOGY 2D Imaging

X-MIND™ UNITY - BOTTOM MOUNTING

XMU BOTT 40 SHORT CONE

HXMU-BS

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 16"

XMU BOTT 80 SHORT CONE

HXMU-BM

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 31"

XMU BOTT 110 SHORT CONE

HXMU-BL

DC Tube Head, Scissor Arm, Timer, 2 Stud Wall Mount, Short Cone, Bottom Mount Extender Arm - 43"

Digital X-ray Systems to integrate into X-MIND™ Unity

CMOS sensor with optic fiber of 1.5 million of pixels (size 1) or 2.21 millions of pixels (size 2), theorical resolution of 25 pl/mm, actual resolution > 12 pl/mm, 3.70m cable. Includes ACE technology (Automatic Control Exposure).

Delivered with a sensor holder, 10 sensor sheaths, and ACTEON imaging Software.

SOPIX® INSIDE SENSOR SIZE 1

S 802 6500

SOPIX (SIZE 1) DIGITAL SENSOR - DIRECT USB2

SOPIX® INSIDE SENSOR SIZE 2 S 802 6501

SOPIX (SIZE 2) DIGITAL SENSOR - DIRECT USB2





INTRAORAL RADIOLOGY Accessories

X-MIND UNITY ACCESSORIES



X-MIND™ Unity short cone W0800018 SHORT CONE: 8" Ø60 (20cm)



X-MIND™ Unity long cone W0800019 LONG CONE: 12" Ø60 (31cm)



X-MIND™ Unity short square cone W0800020 Square SHORT CONE: 8″ Ø60 (20cm)





X-MIND™ Unity remote exposure switch W0800021
REMOTE EXPOSURE SWITCH



WALLPLATE for X-MIND™ Unity installation SSWALLPLATE
Single Stud Wallplate



XMU EXTERNAL LIGHT 115V W0800034

INTRAORAL RADIOLOGY Digital radiology system by phosphor plate



PSPIX^{®2} S 702 0011

Digital imaging plate system

Includes 4 imaging plates (2 x size 1, 2 x size 2), 200 bag & cover for imaging plate (100 x size 1, 100 x size 2), an imaging plate storage box, 1 stylus for touchscreen,1 power supply 100-240 VAC, 1 Ethernet cable, ACTEON imaging software and its license, and a Mac and its license.



STANDARD IMAGING PLATE





	PSPIX ² Set of 2	PSPIX 1 st generation Set of 2
Size 0 (22 x 35mm)	990215	700178
Size 1 (24 x 40mm)	990216	700179
Size 2 (31 x 41mm)	990217	700180
Size 3 (27 x 54mm)	990218	700181

INTRAORAL RADIOLOGY Accessories

IDOT™IMAGING PLATE

Numbered plate for better identification.



	PSPIX ²	PSPIX 1 st generation set of 6
Size 0 (22 x 35mm)		700182
Size 1 (24 x 40mm)		700183
Size 2 (31 x 41mm)		700184
Size 3 (27 x 54mm)		700185

PROTECTIVE BAG FOR IMAGING PLATE



BAG & COVER FOR IMAGING PLATE with integrated cardboard envelope









PSPIX²

Size 0

Size 1

Size 2

Size 3



Size 0



Size 1





PSPIX 1st generation

	PSPIX ² set of 300	PSPIX 1st generation set of 250
Size 0	700546	700338L
Size 1	700547	700339L
Size 2	700548	700340L
Size 3	700549	700341L

PSPIX² ACCESSORIES



Bag & cover for imaging plate start up kit

(set of 100 pieces of size 1 and 100 pieces of size 2)



Imaging plate storage box 990219

(delivered without imaging plates)



Stylus for PSPIX² touchscreen 005132



Cleaning plate for PSPIX² 990243

INTRAORAL RADIOLOGY Accessories



Autoclavable Insertion Guide 608019

Autoclavable insertion guide (Top) for imaging plates

Autoclavable Receptacle 608018

Autoclavable receptacle (Bottom) for imaging plates



Size 4 kit for PSPIX² 990250

1 plastic SUPPORT size 4, 100 protective bags size 4, 3 imaging plates size 3, 1 quick set-up, 1 quick start

Plastic Support Size 4 990251

Plastic case for 3 size 3 plates

BAGS & COVERS S4 - 100 700644

Bag and cover for imaging plates size 4 (set of 100)



PSPIX Positioning Set - Bitewing 700195

Includes five positoning holders with ring for bitewing and one centring devices.



PSPIX Positioning Set - Anterior & Posterior 700196

Includes two positoning holders with ring for anterior teeth, two positoning holders with ring for posterior teeth and two centering devices.



Imaging Plate Storage Box - Gen 1 700194

(delivered without imaging plates)

INTRAORAL RADIOLOGY Digital radiology sensors





CMOS sensor with optic fiber of 1.5 million of pixels (size 1) or 2.21 millions of pixels (size 2), theorical resolution of 25 pl/mm, actual resolution > 18 pl/mm, 3.70m cable. Includes ACE technology (Automatic Control Exposure).

Delivered with a sensor holder, 10 sensor sheaths, and ACTEON imaging software.

SOPIX®2 (SIZE 2) DIGITAL SENSOR S 802 0005

SOPIX² Direct USB2 - comes w/ one sensor holder

$SOPIX^{\otimes 2} \ (SIZE \ 1) \ DIGITAL \ SENSOR$

S 802 0006

 $\overline{\text{SOPIX}^2}$ (SIZE 1) DIGITAL SENSOR - DIRECT USB2 - COMES W/ ONE SENSOR HOLDER

INTRAORAL RADIOLOGY Accessories

ACCESSORIES



Positioning set for SOPIX²/SOPIX series sensors 700042

SOPIX² Positioning Set (5 sets with positioning Ring for: BW, Ant, Post, Endo Ant, Endo Post)



Sopix²/ Sopix sensor holder 926039

Size 1 or Size 2



Sheaths for Sopix²/Sopix

(500 pieces in a single pasteboard) for size 1 SOPIX² series sensors for size 2 SOPIX² series sensors

403026 403027



Sensor Holder Combo Kit ACT-459-1+2

Size 1&2 Combo kit - includes 1 Ring & rod system, 4 size 1 bite blocks, 4 size 2 bite blocks (1 each size Anterior, Posterior, Vertical and Horizontal Bitewings)

Sensor Holder Size 2 Kit ACT-459-2

Size 2 kit - includes 1 Ring & rod system, 8 size 2 bite blocks (2 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

Sensor Holder Size 1 Kit ACT-459-1

Size 1 Kit (includes 1 Ring & rod system, 8 size 1 bite blocks: 2 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

INTRAORAL RADIOLOGY Accessories

Sensor Holder Bite Blocks Size 2 ACT-462

16 pack of bite blocks size 2 (4 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

Sensor Holder Bite Blocks Size 1 ACT-461

16 pack of bite blocks size 1 (4 each: Anterior, Posterior, Vertical and Horizontal Bitewings)

*Sensors, ring and rod not included



3 pack of bite blocks size 1 Anterior *ACT-001*

3 pack of bite blocks size 2 Anterior *ACT-049*

*Sensors, ring and rod not included



3 pack of bite blocks size 1 Posterior *ACT-018*

3 pack of bite blocks size 2 Posterior *ACT-056*

*Sensors, ring and rod not included



3 pack of bite blocks size 1 Horizontal Bitewing *ACT-025*

3 pack of bite blocks size 2 Horizontal Bitewing ACT-063

3 pack of bite blocks size 1 Vertical Bitewing *ACT-032*

3 pack of bite blocks size 2 Vertical Bitewing ACT-070

*Sensors, ring and rod not included

3D imaging system

- Easily upgradeable
- Superior 3D image processing with SHARP and STAR
- Patented Ceph collimation
- Large selection of F.O.V: Ø40x40, 60x60, 80x80 and 110x80

X-MIND™ TRIUM PAN

Control panel on the right Control panel on the left

W0900001 W0900007



OrthoPantomograph (OPG), 115-220 V, With One CMOS Pan 2D sensor for panoramic imaging

- Adult and Child Panoramic programs
- Improved Orthogonality
- TMI section
- Maxillary Sinus exam

ACTEON® Imaging Suite Software, for visualisation and analysis of 2D acquisitions.

Installation and user e-manuals, 10 year warranty.

Not included: Workstation PC and monitor, light.





X-MIND® Trium is a class 3R laser product per IEC 60825-1:2007. Avoid direct eye exposure to laser radiation. Viewing the laser output with magnifying optical instruments (for example, surgical microscope and binocular glasses) may pose an eye hazard and thus the user should not direct the beam into an area where such instruments are likely to be used.

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1

X-MIND™ TRIUM PAN (CEPH READY) UPGRADES

X-MIND Trium Pan Ceph Ready Upgrade Ceph right (no sensor included)
X-MIND Trium Pan Ceph Ready Upgrade Ceph left (no sensor included)

W0900007
W0900007

X-MIND™ TRIUM PAN (3D READY)

Control panel on the right W0900225
Control panel on the left W0900226

UPGRADE 3D

X-MIND Trium Pan 3D Ready Upgrade 3D

W0900227

UPGRADES CEPH

X-MIND Trium Pan 3D Ready Upgrade Ceph right (1 sensor Ceph 2D included)
X-MIND Trium Pan 3D Ready Upgrade Ceph left (1 sensor Ceph 2D included)

W0900018 W0900019









^{*}Specific restrictions apply for 11 x 8cm FOV in Ontario, Canada. Contact ACTEON North America for more information.

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1



X-MIND™ TRIUM PAN CEPH

CEPH left (control panel on the right)
CEPH right (control panel on the left)
CEPH right (control panel on the right)
CEPH left (control panel on the left)

W0900003 W0900009 W0900230 W0900232



OrthoPantomograph (OPG), 115-220V, with one CMOS Ceph 2D sensor used for panoramic imaging and cephalometry.

- Adult and Child Panoramic programs
- Improved Orthogonality
- TMJ section
- · Maxillary Sinus exam
- · Cephalometric exam

ACTEON® Imaging Suite (AIS) Software, for visualisation and analysis of 2D acquisitions.

Installation and user e-manuals, 10 year warranty.

Not included: Workstation PC and monitor, light.

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1





X-MIND™ TRIUM PAN 3D

Control panel on the right Control panel on the left W0900004 W0900010

OrthoPantomograph (OPG), 115-220V, with one CMOS

Ceph 2D sensor used for panoramic imaging and one Flat Panel Detector for CBCT.

- Adult and Child Panoramic programs
- · Improved Orthogonality
- TMJ section
- Maxillary Sinus exam
- CBCT FOV Ø40x40, Ø60x60, Ø80x80, Ø110x80

Workstation PC with standard monitor.

Software preinstalled.

ACTEON® Imaging Suite (AIS) Software, for visualisation and analysis of 2D and 3D acquisitions and treatment planning. Includes five licenses for 3D app, the 3D viewer and implant planning

application.

Installation and user e-manuals, 10 year warranty.

Not included: light.

UPGRADES

X-MIND™ Trium Pan 3D Upgrade Ceph right (1 2D Ceph sensor included)

W0900018

X-MIND™ Trium Pan 3D Upgrade Ceph left (1 2D Ceph sensor included)

W0900019

*Specific restrictions apply for 11 x 8cm FOV in Ontario, Canada. Contact ACTEON North America for more information.

ACTEON IMAGING

EXTRAORAL RADIOLOGY X-MIND® Trium - Upgradeable 3 in 1

X-MIND™ TRIUM PAN CEPH 3D



CEPH left (control panel on the right)
CEPH right (control panel on the left)
CEPH right (control panel on the right)
CEPH left (control panel on the left)

W0900006 W0900012 W0900231 W0900233



OrthoPantomograph (OPG), 115-220V, with one CMOS Ceph 2D sensor used for panoramic imaging and cephalometry and one Flat Panel Detector for CBCT.

- Adult and Child Panoramic programs
- · Improved Orthogonality
- TMJ section
- · Maxillary Sinus exam
- Cephalometry exam
- CBCT FOV Ø40x40, Ø60x60, Ø80x80, Ø110x80

Workstation PC with standard monitor. Software preinstalled.

ACTEON® Imaging Suite (AIS) Software, for visualisation and analysis of 2D and 3D acquisitions and treatment planning. Includes five licenses for 3D app, the 3D viewer and implant planning application.

Installation and user e-manuals, 10 year warranty.

Not included: light.

*Specific restrictions apply for 11 x 8cm FOV in Ontario, Canada. Contact ACTEON North America for more information.

EXTRAORAL RADIOLOGY AIS & Accessories





_

W0911246

ACTEON® IMAGING SUITE

Trium 2D ACQUISITION PC USA

AIS 3D APP DESIGN LICENSE FIRST YEAR	W1000025
3DAPP DESIGN LICENSE ADDL YEAR (REQ W1000025 PRIOR TO PURCHASE)	W1000026
AIS 3D APP STANDARD 1 ELECTRONIC LICENSE	W1000014
AIS 3D APP BASIC 1 ELECTRONIC LICENSE	W1000012
AIS 3D APP BASIC 5 ELECTRONIC LICENSES	W1000013
AIS 3D APP STANDARD 5 ELECTRONIC LICENSES	W1000015
AIS 3D APP PREMIUM 1 ELECTRONIC LICENSE	W1000016
AIS 3D APP PREMIUM 5 ELECTRONIC LICENSES	W1000017
AIS DICOM PACKAGE license (DICOM print license NOT included)	W0900153
DICOM print license	39800751
A CCECCO DIEC	
ACCESSORIES	
Trium POWER CONDITIONER 8A, 120V	61084-03R
Trium 2D ACQUISITION PC AND MONITOR	W0900246

EXTRAORAL RADIOLOGY Accessories

Pan 2D sensor for OrthoPantomograph (OPG) (to add a second 2D sensor or for replacement)	W0900036
Ceph 2D sensor for OrthoPantomograph (OPG) and cephalometry (for replacement)	W0900037
3D sensor Plan captor for 3D OrthoPantomograph (OPG) (for replacement)	W0900038
115V Light Reports X-Rays emission during acquisition	29701225.CW
230V Light Reports X-Rays emission during acquisition	29701226.CW
PC workstation (no monitor) Windows 7, 64 bits. Nvidia Graphic card. Not included: network card, monitor, AIS 3D App	W0900040
PC workstation with network card and dongle key (no monitor) Windows 7, 64 bits. Nvidia Graphic card. 1 Gbit Ethernet Network card. 3D rdongle. Not included: monitor, AIS 3D App.	W0900042 reconstruction
24" PC Monitor	W0900141
Remote trigger & 15m cable	W0900091
Ears protections X-MIND Trium only for arms Ceph	W0900110
XMT CHIN REST, BITE BLOCK, TMJ	W0900080
XMT SLIDING BITE BLOCKS (10pcs)	W0900080
Bite protections X-MIND Trium (500 pcs)	W0900142
X-MIND Trium Self Standing Base	W0900236



X-MIND PRIME 2D PAN

XMP 2D WALL 110-120V: W1200001

COMPLETE SYSTEM KIT: PRIME + WALL FIXING

PLATE: W1200001-08

X-MIND® Prime is a smart space-saving solution, providing high-technology capabilities and simplicity of use. This ingenious unit offers all essential diagnostic tools in a cost effective

package.

X-MIND PRIME 3D CBCT/PAN

XMP 3D WALL 110-120V + WS + MONITOR: W1200005 COMPLETE SYSTEM KIT: PRIME + WALL FIXING PLATE + EXTENDED VIEW + US PLUG: W1200005-09

X-MIND® Prime 3D offers a multitude of acquisition programs thanks to its multiple FOV (5x5cm, 8.5x5cm, 8.5x9.3cm and 12x10cm). Once you select the region of interest of the examination, the settings and dose are automatically adjusted.

X-MIND PRIME COLUMN

Column for X-MIND® Prime for when mounting to the wall is not possible W1200101V

XMP Disposable Bite Protection (100) - 6107110700 XMP Chin Support Kit - 7104011300 XMP TMJ Protective Sleeves (60) - 6107110800

X-MIND PRIME PAN WITH CEPH

X-Mind Prime with Cephalometric Arm W1200009-08

X-MIND PRIME 3D WITH CEPH

COMPLETE SYSTEM KIT: PRIME + WALL FIXING PLATE + EXTENDED VIEW + US PLUG W1200010-09



INTRAORAL CAMERAS Diagnostic tools



Patented SOPRO® innovation using fluorescence and selective chromatic amplification technologies.

SOPROCARE® instantly highlights plaque, gingival inflammation and caries even in their early stages.

- Light (78 g)
- 4 pre-set selections (extraoral, intraoral, CARE, macro)
- Lighting by LED: four white three blue
- Freezing by sensitive touch "SOPROTouch®"



SOPROCARE®

Handpiece *S_950_0002*

Delivered with four SOPROTIPS, a handpiece holder, ACTEON imaging software and its license. Needs to be used with a docking station.



Patented SOPRO® innovation using fluorescence technologies.

SOPROLIFE® Detects caries on occlusal and interproximal surfaces, helping to achieve minimally invasive treatments.

- Light (78 g)
- 4 pre-set selections (extraoral, intraoral, LIFE, macro)
- Lighting by LED: 4 whites (Daylight mode) or 4 blues (diagnostic or traitement mode)
- Freezing by sensitive touch "SOPROTouch®"



SOPROLIFE® Handpiece

5 900 0002

Delivered with four SOPROTIPS, a SOPROLIFE® handpiece holder, ACTEON imaging software and its license. Needs to be used with a docking station.



SOPRO 717 First reveals micro-cracks, micro-lesions, infiltrations – all that human vision can not see alone.

- Light (75 q)
- 3 pre-selections (extraoral, intraoral, macro)
- Lighting by LED (x8)
- Freezing by sensitive touch "SOPROTouch®"



SOPRO® 717 FIRST

Handpiece

S_717_0101

Delivered with a SOPRO® 717 First handpiece holder, ACTEON imaging software and its license. Needs to be used with a docking station.

SOPR **617**

- · Very light (55 g)
- Freezing by sensitive touch "SOPROTouch®"
- No focus adjustment
- Lighting by LED (x8)
- Non-inverted image
- 105° angle of view



SOPRO® 617

Handpiece

5 617 0002

Delivered with a SOPRO® 617 handpiece holder, ACTEON imaging software and its license. Needs to be used with a docking station

L ACTEON

INTRAORAL CAMERAS Docking Stations

DOCKING STATIONS



Mini dock USB 2.0

S_717_1600

Docking station USB2.0

Including one 2.5 meters linking cable and one 0.8 meter repeater cable

Mini Dock U USB 2M Surf6 *S_717_1703*

Mini Dock, 2M cable in Surf, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

Mini Dock U USB 7M Surf6 5 717 1704

Mini Dock, 7M cable in Surf, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

Mini Dock U USB 2M Grey3 S 717 1705

Mini Dock, 2M cable in grey, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

Mini Dock U USB 7M Grey3 5 717 1706

Mini Dock, 7M cable in grey, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately



Linking cable 8ft

922050

With connector on both sides SOPROLIFE®/SOPRO® 717/SOPRO® 617

Linking cable 15ft

922051

With connector on both sides SOPROLIFE®/SOPRO® 717/SOPRO® 617

INTRAORAL CAMERAS Accessories

FOOTSWITCHES



USB Foot Pedal

009009

SHEATHS



Sheaths for SOPROLIFE, 403017-TORIOC20 SOPROCARE or SOPRO 717 First

500 sheaths per box

Sheaths for SOPRO 617

500 sheaths per box

403023

ACCESSORIES



SOPROTIPS for SOPROLIFE® and SOPROCARE®

Protection tip against light. Sterilization by autoclave. (Set of 4 pieces) 990135

ULTRASONICS Ultrasonic generators





NEWTRON® BOOSTER US F60216 CANADA F60206

NEWTRON unit delivered with NEWTRON Handpiece, US power Supply

NEWTRON® BOOSTER WITH INTRO TIP KIT US F60211 CANADA F60201

NEWTRON unit delivered with NEWTRON Handpiece, intro kit (#1, 10Z, H3 tips, three blue torque wrenches, autoclavable plastic box), and US power Supply



NEWTRON® P5 B.LED US F61116 CANADA F61106

NEWTRON P5 B.LED unit delivered with NEWTRON Slim B.LED Handpiece with white LED ring, US power Supply

NEWTRON® P5 B.LED WITH INTRO KIT US F61117 CANADA F61107

NEWTRON P5 B.LED delivered with NEWTRON Slim B.LED Handpiece with white LED ring, intro kit (#1, 10Z, H3 tips, three blue torque wrenches, autoclavable plastic box), and US power Supply

Class IIa medical devices (GMED) - CE 0459

ULTRASONICS Ultrasonic generators





NEWTRON® P5 XS B.LED US F62118 CANADA F62106

NEWTRON P5 XS B.LED unit delivered with NEWTRON Slim B.LED Handpiece with white LED ring, US power Supply

NEWTRON® P5 XS B.LED WITH INTRO KIT US F62117 CANADA F62107

NEWTRON P5 XS B.LED delivered with NEWTRON Slim B.LED Handpiece with white LED ring, intro kit (#1, 10Z, H3 tips, three blue torque wrenches, autoclavable plastic box), and US power Supply

Class IIa medical devices (GMED) - CE 0459

ACCESSORIES



NEWTRON® handpiece

F12281

Autoclavable Handpiece: GII 2-Tone (Grey)
NEWTRON OFM



Nose for NEWTRON® handpiece

E12281



NEWTRON® SLIM B.LED handpiece F12900

NEWTRON Slim B.LED Handpiece with blue LED ring; Compatible with NEWTRON B.LED and NEWTRON XS B.LED devices for use with FLAG fluorescent liquid assistance guide

*F.L.A.G. is currently unavailable in Canada



NEWTRON® SLIM B.LED handpiece

F12905

NEWTRON Slim B.LED Handpiece with white LED ring for NEWTRON B.LED, XS B.LED



NEWTRON® SLIM handpiece without LED

F12901

For NEWTRON® Booster sold before 2016



Nose for NEWTRON® SLIM and NEWTRON® SLIM B.LED handpieces

E12904



NEWTRON® SLIM B.LED blue LED ringBlue LED Ring for NEWTRON Slim B.LED handpiece

(to illuminate FLAG - not available in Canada)



NEWTRON® SLIM B.LED white LED ring

F12915

White LED Ring for NEWTRON Slim B.LED Handpiece

Light Guide

E12903

LED Light Guide for NEWTRON Slim B.LED Handpiece

ULTRASONICS Ultrasonic Accessories



300ml tank

F62004



500ml tank

F62005



F.L.A.G.™ for B.LED

F62209

FLAG 5, 1.5 ML vials. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED handpiece. Not available in Canada.

F.L.A.G.™ for B.LED

F62207

FLAG 100 ML bottle. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED handpiece. Not available in Canada.

F.L.A.G.™ Package

F62206

FLAG Package 100 ML bottle with Blue LED ring. Fluorescent liquid assistance guide adheres to plaque and is illuminated by the blue light of the NEWTRON Slim B.LED.

*F.L.A.G. is currently unavailable in Canada



Large O-Ring for Handpiece *E15019*

Small O-Ring for Handpiece

O-Ring Kit *E15121*

for P Max XS and P5 NEWTRON XS Irrigation Tanks

Medium O-Ring

E15002

E15052

Black Medium O-ring 9 x 1 for NEWTRON Slim and NEWTRON Slim B.LED Handpiece

O-Ring with Tool

F12304

Handpiece 3 gasket O-rings and tool

33 CACTEON EQUIPM

ULTRASONICS Ultrasonic Accessories

Water Filter, Comp	olete	F10386

Replacement water filter cartridge F10389



NEWTRON® LED handpieceAutoclavable LED Handpiece for P Max XS, P5 LED, P5 XS LED, SP NEWTRON LED and Piezotome 2 / Implant Center 2



Light RingFor NEWTRON® LED handpiece

F12605

Light Guide E12674Light guide for the NEWTRON® LED Handpiece



Blue Plastic Color Coded Wrench F81302



NEWTRON® wrenchBlue color coded torque wrench

F81330



Universal wrench

F00406



Autoclavable plastic box (round) F12351
Clear Cover GII (Empty) Autoclavable Box - Holds
four tips total - Requires round torque wrenches
to hold the tips (not included) - F81330 and any

previous round torque wrench versions

Autoclavable plastic box


F12350


Clear Cover GII (Empty) Autoclavable Box - Holds four tips total (not included)



Autoclavable box F00143

Metal Blue Autoclavable Box - Holds 6 tips (not included)











IP1 F02121 **Implant** abutments and wide threads



F02122 threads, left-angled



F02123 Medium implant Medium implant threads, right-angled



IP3L F02124 Narrow implant threads, left-angled



F02125 Narrow implant threads, right-angled

perio**m**aintenance



F01001 Short probe



F01004 Long probe



F02162 Premolars & molars. left-angled



TK2-1R F02161 Premolars & molars. right-angled

perioSoft



PH1 (x4) F00702 Anterior teeth



PH2L (x4) F00705 Premolars & molars. left-angled



PH2R (x4) F00706 Premolars & molars. right-angled

erio Precision



P₂L F00090 Narrow pockets debridement, left-angled



P2R F00091 Narrow pockets debridement. right-angled

endodontics



F43710

Ø 0.1/length 21mm





F43712

Ø 0.1/length 25mm



K15-21 (x4) F43715 Ø 0.15/length 21mm



K15-25 (x4) F43717 Ø 0.15/length 25mm



K25-21 (x4) F43725 Ø 0.25/length 21mm



K25-25 (x4) F43727 Ø 0.25/length 25mm



K30-21 (x4) F43730 Ø 0.3/length 21mm



K30-25 (x4) F43732 Ø 0.3/length 25mm

irrisafe





IRR 20-25 (x4) F43808 Ø 0.2mm length 25mm

IRR 25-21 (x4) F43805 Ø 0.25mm length 21mm

IRR 25-25 (x4) F43806 Ø 0.25mm length 25mm

endosurgery



P14D F00106 Universel



P15LD F00107 Left-angled




P15RD F00108 Rightangled




S12-70D F00118 70° Universal

















Class IIa medical devices (GMED) - CE 0459







scaling

Scaling kit

F00900

Tips - #1, #2, 10X, 10Z with 4 torque wrenches and Autoclavable Box



NEWTRON

NEWTRON Intro Kit

F00932

Tips - #1, 10Z, H3, three torque wrenches and Autoclavable Box



perio Precision

Perio Precision kit

F00939

Tips - TK1- 1S, P2L, P2R, 3 toque wrenches and Autoclavable Box



periodontics

Perio kit

F00902

Tips - H3, H4L, H4R with 3 torque wrenches and Autoclavable Box



perio*m*aintenance

Bio-Film Disruption Kit

F00904

Tips - TK1- 1S /TK1- 1L / TK2- 1R / TK2- 1L with 4 torque wrenches and Autoclavable Box



implant Protect

Implant Protect kit

F02120

Tips - IP1, IP2L, IP2R, IP3L, IP3R with universal wrench and Autoclavable stand



periodontics

Perio Soft Kit

F00906

Tips - 4 pk PH1, 4 pk PH2L, 4 pk PH2R with 3 blue torque wrenches and Autoclavable Box. Only available until supplies last.



excavus

Excavus® kit

F00739

Tips - EX1, EX2, EX3, EXR, EXL with 1 Universal Wrench, and Autoclayable Stand



endosuccess Retreatment

Endo Success™ Retreatment kit

F00737

ET18D, ETBD, ET20, ET25S, ET25 and ETPR with 1 Universal Wrench and Autoclavable Stand



endosuccess Canal Access Prep

Endo Success™ Canal Access Prep kit

F88180



CAP1, CAP2, CAP3 with universal wrench and Autoclavable stand



endosuccess Apical Surgery

Endo Success™ Apical Surgery kit

F00069

AS3D, AS6D, AS9D, ASLD, ASRD Tips with 1 Universal Wrench and Autoclavable Stand

ULTRASONICS Tip kits





Perfect Margin® Rounded kit

F00738

Tips - PM1, PM2, PM3, PM4 with 1 Universal Wrench and Autoclavable Stand





Perfect Margin® Shoulder kit

F00736

Tips - PMS1, PMS2, PMS3, PM4 with 1 Universal Wrench and Autoclavable Stand



perfectmargin Veneers



Perfect Margin® Veneers kit

F02020

Tips - PMV1, PMV2, PMV3, PMV4, PMV5, PMV6 with universal wrench and Autoclavable stand



Micro Retro Kit

F00912

Tips - P14D, P15LD, P15RD with 3 Wrenches and Autoclavable Box



AIR-N-GO® easy

The new generation of dual purpose Air Polisher: SUPRA and PERIO

AIR-N-GO® easy with KaVo* adapter	F10130
AIR-N-GO® easy with Sirona* adapter	F10131
AIR-N-GO® easy with W&H* adapter	F10132
AIR-N-GO® easy with Bien-Air* adapter	F10133
AIR-N-GO® easy with Midwest* adapter	F10135

Delivered with a convertible AIR-N-GO® easy handpiece, a SUPRA nozzle, 1 intro kit of 20g bags of CLASSIC powders, 1 AIR-N-GO® easy maintenance kit (replacement o-rings, syringe of silicone grease with brushes, cannula, cleaning drill).

Class IIa medical devices (GMED) - CE 0459





 $^{{\}rm *KaVo},$ Sirona, W&H, Bien-Air and Midwest are registered trademarks.

AIR-N-GO® EASY ACCESSORIES

FOR AIR-N-GO EASY ORIGINAL PURCHASED PRIOR TO 2015



AIR-N-GO® easy PERIO KIT

F10120

Turns your handpiece into a Perio system to treatperiodontal disease (tooth & implants) on pockets from 3 to 8mm.

Delivered with a AIR-N-GO easy nozzle and a bottle of AIR-N-GO PERIO glycine powder.

Class IIa medical devices (GMED) - CE 0459



AIR-N-GO® SUPRA nozzle

F10125

For all prophylaxis supra-gingival treatments.



AIR-N-GO® PERIO easy nozzle

F10126

Turns your handpiece into a Perio system to treat non surgical periodontal disease (tooth & implants) in pockets from 3 to 8mm.



AIR-N-GO® PERIO nozzle

F10127

Turns your handpiece into a Perio system to treat non surgical and surgical periodontal disease (tooth & implants) in pockets from 8 to 10mm.



AIR-N-GO® PERIO Maintenance nozzle

F10128

Turns your handpiece into a system adapted for periodontal maintenance in pockets until 4mm.



Air-N-Go Easy O-ring kit

F10121

Air-N-Go Easy o-ring kit: includes silicone grease syringe, (3) brushes, o-rings

POWDER - BOTTLES

Class IIa medical devices (GMED) - CE 0459



AIR-N-GO® CLASSIC powder

4 bottles (250g)

F10250
F10251
F10252
F10253
F10254



4 flavors AIR-N-GO® CLASSIC powder

4 bottles (250g)

F10255



AIR-N-GO® PERIO powder

F10270

3 bottles (160g)

POWDER - STICKS

Class IIa medical devices (GMED) - CE 0459



AIR-N-GO® CLASSIC powder

50 sticks (20g)

Peppermint	F10256
Lemon	F10257
Raspberry	F10258
Cola	F10259
Neutral	F10260

AIR-N-GO® ACCESSORIES

	Ø0.5 mm AIR-N-GO® SUPRA probe	F10207
	Ø1 mm AIR-N-GO® SUPRA probe	F10208
	AIR-N-GO® PERIO probe	F10206
•0	AIR-N-GO® o-rings kit	F10205
	AIR-N-GO® Easy Front Body Autoclavable front body for the AIR-N-GO E	E10207 Easy



OTHER ACCESSORIES



Blue ProphyPenDisconnectable polisher handpiece (bayonet) delivered with a 135° titanium nozzle and an autoclavable universal wrench.



White ProphyPen F10557
Disconnectable polisher handpiece
(bayonet) delivered with a 135° titanium nozzle and an autoclavable universal wrench.



Titanium nozzle 135° Ø0.7 mm

GRACEY CURETTES



GRACEY 1-2 Ref. 402.01SI GRACEY 13-14 Ref. 402.13SI GRACEY 5-6 Ref. 402.05SI GRACEY 15-16 Ref. 402.15SI GRACEY 11-12 Ref. 402.11SI GRACEY 17-18 Ref. 402.17SI

UNIVERSAL CURETTES



 BARNHART 5-6
 Ref. 408.05SI
 COLUMBIA 4L-4R
 Ref. 405.04SI

 LANGER 5-6
 Ref. 406.06SI
 COLUMBIA 13-14
 Ref. 405.13SI

SICKLE SCALERS



M23 Ref. 404.23SI H6-H7 Ref. 406.07SI

COMBOS



H5-L5 Ref. 406.05SI N128-L5 Ref. 406.08SI N128-L5 mini Ref. 406.09SI

DIAGNOSTICS - EXPLORERS



RESTORATION - EXCAVATORS & BURNISHERS



RESTORATION - COMPOSITE & CEMENT SPATULAS



PERIODONTICS - EXPROS



ENDODONTICS - ROOT CANAL PLUGGERS



CASETTES



Easyclip Cassette 8 Instruments Dimensions 180x136x30mm Material: PSU Handy and light (90g) Sterilizable up to 180°C Instruments not included



Sterilization Cassette 16 Instruments Dimensions 285x92x37mm Material: stainless steel Yellow silicone + 2 metal clips Sterilizable up to 180°C Instruments not included

PROPHYLAXIS KIT



Prophylaxis Kit

400CSS.SI

Bliss Prophy kit includes: 1 M23, 1 Gracey 11-12, 1 Gracey 13-14, 1 H6-H7, 1 Columbia 13-14, 1 Bliss Mirror Handle, 1 #5 Pure Reflect Mirror, 1 Probe, 1 Autoclavable cassette

233.20

233.30







Pure Reflect Micro Mirror #5 (CSS) Box of 12 mirrors, single sided 9x3 mm

MP4624PH



Pure Reflect Micro Mirror #4 (CSS) Box of 12 mirrors, single sided 7x2 mm

MP4622PH



Ultrasonic Bone Surgery provides Less Trauma and Less Drama for:

- Atraumatic Extractions
- Bone grafting
- · Sinus elevations
- Ridge Splitting
- Accelerated Orthodontics



PIEZOTOME® CUBE

US F50108 CANADA F50117

Piezotome CUBE with power cord, Piezotome CUBE handpiece, torque wrench, foot pedal, irrigation bracket, and five (5) sterile disposable irrigation lines


Class IIa medical devices (GMED) - CE 0459

NOTE

Some accessories adapted to Piezotome® 2, Piezotome® Solo, Piezotome® CUBE and ImplantCenter™ 2 are not compatible with the first generation Piezotome® and ImplantCenter™, and vice-versa

- References compatible with 1st generation devices are indicated by a blue square :
- References compatible with 2^{nd} and 3^{rd} generation devices are indicated by a grey square \Box

PIEZOTOME® CUBE ACCCESSORIES



CUBE

IMPLANT CENTER 2



FUNCTIONS

- I-Surge™: implantology
- Piezotome[®]: bone surgery
- NEWTRON®: conventional treatments (perio, endo, ...)

IMPLANT CENTER 2 LED: STARTER KIT

F27168

Delivered with i-Surge LED motor, i-Surge LED motor cord, power cord, a progressive footswitch, 2 brackets, 2 handpiece holders and 2 sterilizable irrigation lines, 30 perforators, an Extraction II pack *Contra Angles not supplied

This reference does not include the PIEZOTOME 2 handpiece or tips.

IMPLANT CENTER 2 LED STARTER KIT WITH LED SURGICAL F27168-CA CONTRA ANGLE (20:1)

Delivered with i-Surge LED motor, i-Surge LED motor cord, power cord, a progressive footswitch, 2 brackets, 2 handpiece holders and 2 sterilizable irrigation lines, 30 perforators, an Extraction II pack and LED Surgical Contra Angle (20:1) - Bien Air

This reference does not include the PIEZOTOME 2 handpiece or tips.

- References compatible with 1st generation devices are indicated by a blue square :
- References compatible with 2nd and 3rd generation devices are indicated by a grey square

ACCESSORIES FOR IMPLANT CENTER™



Implant Center 2 Foot Control F27170

Implant Center 2 Progressive Foot Pedal



ImplantCenter™ Bracket

E27113

Irrigation for use with Implant Center 2, Piezotome Solo, Piezotome 2





I-Surge™ LED motor *F27210* I-SURGE LED MOTOR for use with the

Implant Center 2 LED only *Contra Angles not supplied



i-Surge™ LED Motor Cord

F27211

Cord for use with the i-Surge LED Motor on Implant Center 2 LED only *Contra Angles not supplied

SURGERY Accessories

IMPORTANT

The handpieces are equipped with keyway systems to prevent a Piezotome® insert from being screwed onto a NEWTRON® handpiece or a classic insert (scaling, periodontics, endodontics, etc.) from being screwed onto a Piezotome® handpiece.



NEWTRON® Handpiece

F12281

Autoclavable Handpiece: GII 2-Tone (Grey) NEWTRON OEM. Compatible with Piezotome 1 and Implant Center 1. Need F57305 to connect.



Piezotome® LED handpiece + cord F12815 Compatible with Piezotome Solo,

Piezotome 2 and Implant Center 2 device



F12821 **LED Ring** Light Ring for the Piezotome 2 LED Handpiece



NEWTRON® LED handpiece

F12609

Autoclavable LED Handpiece for Piezotome 2 / Implant Center 2. Need F12830 to connect.



NEWTRON® LED cord

F12830

Sterilizable and disconnectable cord to connect the NEWTRON LED handpiece to the Implant Center 2 or Piezotome 2

• References compatible with 1st generation devices are indicated by a blue square :



• References compatible with 2nd and 3rd generation devices are indicated by a grey square 🔳



Sterilization Box

F27156

Piezotome (PZI, PZII) Sterilization box: can hold one handpiece with cord, 2 Tip Kits and 1 Wrench (items are not included)



Irrigation Kit (Sterile Disposable)

Kit includes 10 sterile disposable lines and cartridge.



Irrigation Kit (Autoclavable)

F57374

Kit includes 1 autoclavable irrigation line and 10 sterile disposable perforators. (For use with all I-Surge, Piezotome / Piezotome 2, Implant Center / **Implant Center 2)**



Perforators (Disposable)

F57372

10 Disposable perforators for use with the F57374 Irrigation lines







Irrigation Clip

E57353

Single Irrigation Clip for Surgical Ultrasonic and Micro Motor Devices





Piezotome 2 Foot Control

F57404

Piezotome 2 Progressive Foot Pedal





Piezotome SOLO: Foot Pedal

F57504



Bracket

E27113

Irrigation for use with Implant Center 2, Piezotome Solo, Piezotome 2



Piezotome LED Handpiece

F4204F

Piezotome 2, Piezotome Solo LED and Implant Center 2 LED Handpiece and Cord: Nondisconnectable assembly



Piezotome 2 Nose Cone

E12813

Piezotome 2 stainless steel nose cone

.....

LED Light Guide

E12812

Light Guide fo the Piezotome 2 LED Handpiece



- References compatible with 1st generation devices are indicated by a blue square
- \bullet References compatible with 2^{nd} and 3^{rd} generation devices are indicated by a grey square \blacksquare



Sinus Lift-2 Kit

Tips - SL1-2, SL2-2, SL3-2, SL4-2, SL5-2, Autoclavable Stand (These tips will only work with the Piezotome HP) for lateral approach sinus elevation.

Sinus Lift-2 Kit



Sinus Lift-2 tips

SL1-2	F87511
SL2-2	F87512
SL3-2	F87513
SL4-2	F87514
SL5-2	F87515



Intralift[™]-2 Kit

Tips - TKW1, TKW2, TKW3, TKW4, TKW5, Autoclavable Stand to perform internal sinus elevation.

Intralift[™] 2 Kit **F87536**



Intralift™-2 tips

TKW1-2: Conical tip ø1.35mm (pilot drilling) TKW2-2: Cylindrical tip ø2.1mm

(preliminary drilling)

TKW3-2: Cylindrical tip ø2.35mm (preliminary drilling)

TKW4-2: Cylindrical tip ø2.8mm (secondary drilling)

F87531

F87532

F87533

F87534



TKW5-2 Non cutting "trumpet" tip (sinus membrane elevation)



EXTRACTION 2 PACK

Sterilization Box, Piezotome 2 Handpiece, Piezotome 2 Cord, Extraction 2 Kit and Easytorque Wrench for atraumatic extractions

Extraction 2 Pack F87540



Extraction-2 Kit

Tips - LC1-2, LC2-2 (2), LC2R-2, LC2L-2, Ninja-2, Autoclavable Stand and Universal Wrench for atraumatic extractions.

Extraction-2 Kit **F87546**



Extraction-2 tips

LC1-2 LC2-2 LC2L-2 LC2R-2 F87507 F87542 F87543 F87544



Ninja®-2



Crown Extension-2 Kit

BS6-2, CE1-2, CE2-2, CE3-2 Tips with 1 Universal Wrench and Autoclavable Stand to perform crown lengthening procedures.

CE-2 Kit F87554



Crown Extension-2 tips

CE1-2 - Diamond-coated spherical tip ø1.75mm

CE2-2 - Diamond-coated spherical tip ø1.20mm

CE3-2 - Diamond-coated straight tip on 5mm (laser mark placed at 3mm: biologic width) F87551

F87553



BS6-2

F87506



Specialist Upgrade Tip Kits

Bone Surgery-2 Tip Kit, Sinus Lift-2 Tip Kit, IntraLift-2 Tip Kit, Clinical Book and DVD



Crest Splitting-2 Kit

CS1, CS2, CS3, CS4, CS5, CS6, tips an autoclavable metal support, an autoclavable universal wrench. CS-2 Kit **F87567**



Crest Splitting tips

CS1-2 CS2-2 CS3-2 CS4-2 F87561 F87562 F87563 F87564



CS5-2 CS6-2 F87565 F87566



Piezocision™ Kit

Delivered with PZ1, PZ2L, PZ3R, PZ3 tips, an autoclavable metal support and an autoclavable universal wrench.

F87576

F87571

F87572



Piezocision™ tips

PZ1 - Anterior corticotomies.

PZ2L (left oriented)
Lateral corticotomies on patient's right vestibular arch
To be used with pendulum motion

PZ2R (right oriented)
Lateral corticotomies on patient's left vestibular arch
To be used with pendulum motion

PZ3 - Anterior corticotomies for root proximity concerns

F87574





Essential-2 Kit

BS1S II, BS4 II, LC2 II, SL1 II, SL2 II, SL3 II tips, an autoclavable metal support, an autoclavable universal wrench

F87528

NEWTRON® RETRO PACK

Delivered with a sterilization box, an autoclavable NEWTRON® handpiece, an autoclavable NEWTRON® cord, a NEWTRON® Retro kit and a Piezotome® autoclavable wrench.

- References compatible with 1st generation devices are indicated by a blue square
- References compatible with 2^{nd} and 3^{rd} generation devices are indicated by a grey square

IMPORTANT

NEWTRON® kits are compatible with the NEWTRON® and NEWTRON® LED handpieces.



NEWTRON® Perio Kit

Tips - H1, H2L, H2R, H3, H4L, H4R, Autoclavable Stand





NEWTRON® Implant Protect Kit

Tips IP1, IP2L, IP2R, IP3L and IP3R, an autoclavable metal support, an autoclavable universal wrench.





NEWTRON® Micro-Retro Kit

Tips - P14D, P15LD, P15RD with 3 Yellow Wrenches and Autoclavable Box





NEWTRON® Retro Kit

Tips - S12-70D, AS6D, ASLD, ASRD, Autoclavable Stand

F87326



Autoclavable Stand - NEWTRON Tips

NEWTRON- Silver Metal Autoclavable Stand - Holds 6 tips (not included)

N88449



MINILED SUPERCHARGED

Mini LED™ SuperCharged

F05214

2,200 mWcm2 Output, (Modes: Full, Pulse & Ramp) Handpiece, Battery, Charger, Power Supply, Light Guide 7.5mm, Oval Light Shield (Cordless)

Mini LED™ SuperCharged handpiece

F05217





Mini LED™ Ortho 2

F05220

3,000 mWcm2 Ouput (Modes 4, 8, 16 & 32 seconds) Handpiece, Battery, Charger, Power Supply, Light Guide 5.5mm, Oval Light Shield (Cordless)



Mini LED™ Black

F05400

Mini LED Black handpiece, 7.5mm light guide and light shield. Handpiece requires 3 AA batteries (Batteries are included)


POLYMERIZATION Curing lights





Mini LED™ Standard

F02611



Wavelength: 420 - 480nm

Li-Ion battery

ruur

Handpiece: anodized aluminium

Mini LED™ Handpiece

F02530

Mini LED Handpiece for Battery operated unit

POLYMERIZATION Accessories

For Mini LED™ SuperCharged, Mini LED™ Ortho 2

	Mini LED Battery Charger Base	F02510
(Mini LED Battery Lithium Ion	F02520
7-A 7-6	Power supply with US plug for Mini LED	F02544
	Light Shield Light shield for Mini LED, SuperCharged, Or Black	F05407 tho2 &
	Oval Light Shield Oval Light Shield for Standard Mini LED	F02555

For Mini LED™, Mini LED™ SuperCharged, Mini LED™ Ortho 2



PROSTHESIS & AESTHETICS Gingival Retraction Paste





EXPASYL™ EXACT

Effective and atraumatic gingival retraction for traditional and digital impressions, Class II and V Restoration, sealing/bonding, and Implantology. EXPASYL Exact Capsules work with virtually any applicator gun, and each capsule contains one dose per patient.



EXPASYL[™] EXACT APPLICATOR GUNApplicator gun for Expasyl Exact capsules

261036



EXPASYL™ EXACT CAPSULE - 20 CAPSULE PACK 261010

EXPASYL™ EXACT CAPSULE - 50 CAPSULE PACK 261011

*EXPASYL™ Exact is currently unavailable in Canada. See next page for standard EXPASYL™

Manufacturer: Produits Dentaires Pierre Rolland - FRANCE

PROSTHESIS, PROPHYLAXIS & PERIODONTICS Gingival Retraction, Hemostasis and Splinting





EXPASYL™ INTRO KIT

Expasyl Intro Kit (Applicator Gun, 6 Capsules, 12 Applicator Tips) System for temporary gingival retraction.

294100



1 MANUAL APPLICATOR

Expasyl Application Gun

260900



EXPASYL[™] CAPSULE - 20 PACK (STANDARD)

261030

261001

EXPASYL™ CAPSULE - 20 PACK (STRAWBERRY)

EXPASYL!

EXPASYL™ APPLICATOR TIP - BOX OF 40

261040

EXPASYL™ APPLICATOR TIP - BOX OF 100 261005

HEMOSTASYL

Hemostatic dressing with mechanical action

HEMOSTASYL™ INTRO KIT - 2 Syringes, 40 Tips

HEMOSTASYL™ SYRINGE - 2 Syringes

HEMOSTASYL™ APPLICATOR TIP - Box of 40 Tips

246200

246050

246100

*HEMOSTASYL™ is currently unavailable in Canada



Splinting ribbons

BIOSPLINT SPLINTING RIBBONS - REFILL

201661

BIOSPLINT Ribbons: 8 Splinting Ribbons

*BIOSPLINT is currently unavailable in Canada

Manufacturer: Produits Dentaires Pierre Rolland - FRANCE









Disposable tips for air/water syringes



Box of 250 non sterile tips



White	20166
Blue	20166
Yellow	20166
Green	20166
Pink	20166

Box of 250 non sterile tips



White/Mint 201700



RISKONTROL ADAPTERS

Adapter ADEC II - fits Riskontrol Perfect only	200195
Adapter BELMONT	200158
Adapter Dental EZ	200164
Adapter DCI - fits Riskontrol Perfect only	200152
Adapter Siemens 4000	200178
Adapter Sirona	200154
Adapter KAVO	200171
Adapter Luzzani - fits Riskontrol Perfect only	200150

Manufacturer: Produits Dentaires Pierre Rolland - FRANCE All information for a perfect use of these medical devices are indicated on the website: www.acteongroup.com. Products reserved for professional dental use.



MODULES



SP NEWTRON® 4K5V MODULE with connectors (10 & 4 pin)

F55297



SP NEWTRON® 4K3, 7V MODULE

F55287

Module only. Does not contain harness adapter (available seperately if needed).

PIEZO ULTRASONIC MODULE F55480

SP NEWTRON universal w/ connector, AHP. handpiece tubing, Potentiometer. Note: Does not include Air-Electric switch or power supply. Items are available separately.

PIEZO ULTRASONIC MODULE WITH **HYGIENE TIPS** F55331

Reference F55480 with Autoclavable Box, 3 color coded torque wrenches, Tips: #1, 10Z, H3. Note: Does not include Air-Electric switch or power supply. Items are available separately.





MODULES

PIEZO ULTRASONIC MODULE LED F55481

SP NEWTRON LED (Universal), w/ LED handpiece, HP tubing Potentiometer. Note: Does not include Air-Electric switch or power supply. Items are available separately.

PIEZO ULTRASONIC MODULE LED W/ HYGIENE TIPS F55521

Reference F55481 with Autoclavable Box, three color coded torque wrenches, Tips: # 1, 10Z, H3. Note: Does not include Air-Electric switch or power supply. Items are available separately.

HANDPIECES



NEWTRON® handpiece autoclavable, light grey/dark grey

F12281



NEWTRON[®] **LED** handpiece autoclavable

F12609



NEWTRON® SLIM B.LED handpiece with white LED ring

F12905

HANDPIECES ACCESSORIES



ACTEON INTEGRATED _

ULTRASONICS Accessories

ACCESSORIES





Potentiometer F00524
Potentiometer with wire harness F55213
Potentiometer Complete F55213K

Potentiometer, wire harness, knob & sticker

Potentiometer with wire harness* F55213PC1

- Pelton & Crane

Potentiometer with wire harness* F55213PC2

- Pelton & Crane

*Tubing reference depends on Pelton delivery head series

Knob and Sticker for the SP NEWTRON

Knob and Sticker for the SP NEWTRON

TRANSFORMER 36V 2.78A+CORD DTHPU101-36C13M

Power supply for SP NEWTRON + Mini LED + BUILT-IN SOPRO Dock Needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

Air Electric Switch

E00596-OEM

For automatic operation needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

Water Adapter

M00100

Water Adapter for SP NEWTRON

Water Solenoid Kit

M00103

Water Solenoid Kit for installation without water block

Air Bleed Pack

M00101

Air Bleed Pack for SP NEWTRON

ACTEON INTEGRATED

F00526

Handpiece Tubing A-dec (61")	F00184
Handpiece Tubing A-dec GREY3 (61")	F01153
Handpiece Tubing A-dec (81")	F00188
Handpiece Tubing A-dec GREY3 (81")	F01154
Handpiece Tubing LED Surf6 (61")	F00616
Handpiece Tubing LED Surf6 (81")	F00629
Handpiece Tubing LED Grey3 (61")	F00632
Handpiece Tubing LED Grey3 (81")	F00633
Handpiece Tubing - Light grey silicone	F00195
Handpiece Tubing - Light grey silicone with connectors	F00197

Handpiece Tubing - Light grey silicone *F00197PC1* with connectors - Pelton and Crane. Tubing reference depends on Pelton delivery head series.

Handpiece Tubing - Light grey silicone *F00197PC2* with connectors - Pelton and Crane. Tubing reference depends on Pelton delivery head series.

Handpiece Tubing Slim LED G3 61"	F01205
Handpiece Tubing Slim LED G3 81"	F01206

Handpiece Tubing- Light grey silicone F00631 for use with SP NEWTRON LED



Mini LED™ OEM (Built In)

Mini LED Module, Light Grey Cord, Handpiece, Light Guide 7.5mm, Protective Light Shield (Fast, Pulse & Ramp) *Note - Requires a 24 VAC Transformer.

with light grey cord

F02720

Mini LED™ OEM handpiece

F02710

O ACTEON INTEGRATED





Mini LED™ SuperCharged OEM (Built In)

2,200 mWcm2 Output, (Modes: Full, Pulse & Ramp) Handpiece, Module, Light Grey Cord, Light Guide 7.5mm, Oval Light Shield) Note - Requires a 24 VAC Transformer.

with light grey cord

F05205

Mini LED™ SuperCharged OEM handpiece

F05218



	OEM module	F02730
	TRANSFORMER 36V 2.78A+CORD Power supply for SP NEWTRON + Mini LED + E SOPRO Dock Needed for retro fitted NEWTROI or CAMERA OEM installation	BUILT-IN
	Air Electric Switch For automatic operation needed for retro fitte MINI LED or CAMERA OEM installation	E00596-OEM d NEWTRON,
0	Mini LED Handpiece Tubing OEM - light grey Mini LED Handpiece Tubing OEM - light grey - Pelton & Crane	F02721 F02721PC1
	ADEC [*] Surf 6 cord	F02326
	ADEC* 500 Gray 3 (medium grey) cord	F02327

*Adec is a registred trademark

INTRAORAL CAMERAS Docking Stations

BUILT-IN DOCKING STATIONS

Mini Dock U USB 2M Surf6 S 717 1703

Mini Dock, 2M cable in Surf, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

Mini Dock U USB 7M Surf6 S 717 1704

Mini Dock, 7M cable in Surf, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

Mini Dock U USB 2M Grey3 5 717 1705

Mini Dock, 2M cable in grey, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

Mini Dock U USB 7M Grey3 S 717 1706

Mini Dock, 7M cable in grey, USB cable, and ACTEON imaging Software. Air-electric switch not included - must be ordered separately

LINKING CABLES FOR BUILT-IN DOCKING

White Linking cable 21ft

922067

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Grey Linking cable 21ft

922070

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073PC1

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - cable needed depends on Pelton delivery head series number.

Light Grey Linking cable 21ft

922073PC2

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - cable needed depends on Pelton delivery head series number.



INTRAORAL CAMERAS Linking cables



White Linking cable 21ft

922067

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Grey Linking cable 21ft

922070

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073

FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit

Light Grey Linking cable 21ft

922073PC1

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - linking cable needed depends on Pelton delivery head.

Light Grey Linking cable 21ft

922073PC2

Pelton & Crane FOR BUILT-IN DOCK

SOPRO camera linking cable to build into the arm of the dental unit - linking cable needed depends on Pelton delivery head.

SURF6 LINKING CABLE

922118

FOR DOCK MU USB w/swivel connector (7m/23ft) w/clip

SURF6 LINKING CABLE

922117

FOR DOCK MU USB w/swivel connector (2m/6ft) w/clip

GREY3 LINKING CABLE

922127

FOR DOCK MU USB w/swivel connector (7m/23ft) w/clip

GREY3 LINKING CABLE

922126

FOR DOCK MU USB w/swivel connector (2m/6ft) w/clip

SURF6 LINKING CABLE

922139

FOR DOCK U USB w/swivel connector (7m/23ft) w/clip

SURF6 LINKING CABLE

22122

FOR DOCK U USB w/swivel connector (2m/6ft) w/clip

GREY3 LINKING CABLE

FOR DOCK U USB w/swivel connector (7m/23ft) w/clip

GREY3 LINKING CABLE

922140

FOR DOCK U USB w/swivel connector (2m/6ft) w/clip

INTRAORAL CAMERAS Linking cables



GREY LINKING CABLE 21ft

922163

FOR BUILT-IN U USB DOCK SOPRO camera linking cable to build into the arm of the dental unit

GREY3 LINKING CABLE

922141

FOR DOCK U USB w/swivel connector (7m/23ft) w/clip

LIGHT GREY LINKING CABLE 21ft

922143

FOR BUILT-IN U USB DOCK

SOPRO camera linking cable to build into the arm of the dental unit

8M USB Cable

123540

USB CABLE 26 ft for Dock MUSB reference

TRANSFORMER 36V 2.78A+CORD


DTHPU101-36C13M

Power supply for SP NEWTRON + Mini LED + BUILT-IN SOPRO Dock Needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

Air Electric Switch E00596-OEM

For automatic operation needed for retro fitted NEWTRON, MINI LED or CAMERA OEM installation

A GITFON INTEGRATED



RELAX. WE'VE GOT YOU COVERED

ACTEON® offers a One-Year, No-Questions-Asked Warranty on all its devices.* If anything goes wrong within the first year, we will fix it for free.**

MINIMALLY
INVASIVE
SOLUTIONS

^{*} Excludes all consumable products including but not limited to: ultrasonic tips, wrenches, imaging plates, prophy powder, air-water syringes, curing light guides.
** ACTEON reserves the right to change this policy without notice or modify accordingly due to

^{**} ACTEON reserves the right to change this policy without notice or modify accordingly due to malicious or deliberate damage.