

**D.LUX+**

Dual-Wavelength Cordless Curing Light

**BUY 1, GET 1 OR
20 DIAFIL/40 DIAFIL FLOW FREE***

Scan for video



- 1 Powerful 2,400 mW/cm² output
- 2 Dual wavelength: 385 - 515 nm
- 3 360° swivel light probe
- 4 5 versatile curing modes
- 5 Built-in Light Intensity Sensor
- 6 High capacity 2,600 mAH battery
- 7 Wireless inductive charger

**[4008-1001] D-Lux+ Regular Kit**

Contains: D-Lux+ hp & charger, light probe and protector, disposable sheaths (200/box), power supply (adapter)

*DiaFil syringe refills or DiaFil Flow syringe refills. Choices of shades.

DIAFIL™

Light-Cured Nanohybrid Composite

**BUY 2, GET 1 FREE OR
BUY 4, GET 4 FREE!**

*of the same value

- 1 High tensile strength & low shrinkage
- 2 Easy to place/sculpt and does not slump
- 3 Does not stick to instruments
- 4 Adapts well to margins
- 5 Wear & high fracture resistant
- 6 12 shades with remarkable color stability
- 7 Durable restoration w/ long-lasting gloss shine

**SYRINGE REFILL (1 x 4g)**

Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20, A30

**CAPSULE REFILL (20 x 0.25g)**

Shades: A1, A2, A3, A3.5, B2

DIA-X BOND™**Universal**

One-Bottle Universal Bonding Agent

BUY 1, GET 1 FREE!

Scan for peer review



- 1 Versatile with all etching techniques
- 2 All direct and indirect restorations
- 3 Light-, self- & dual-cure compatible
- 4 Virtually no post-op sensitivity
- 5 Advanced 10-MDP monomer
- 6 Convenient: true 1-bottle bond
- 7 Uncompromising bond strength
- 8 Stable at room temperature

**[A2001-2201] Dia-X Bond Universal**

Contains: 1 x 5 ml bottle

DIAFIL™ FLOW

Light-Cured Flowable Composite

**BUY 2, GET 1 FREE OR
BUY 4, GET 4 FREE!**

- Excellent color adaptability
- Pit & fissure sealant
- Excellent tensile strength and low shrinkage
- Highly radiopaque

**REFILL PACK (2g syringe & 10 tips)**

Shades: A1, A2, A3, A3.5, B2, B3, C2, A20, A30.

DIAETCH™

37% Phosphoric Acid Etching Gel

BUY 2, GET 1 FREE!

- Fast! 15-second etch
- Precise! 25 gauge pre-bent tips
- Washes off easily with no residue
- Beautiful blue color
- Ideal consistency and biocompatible

**[A2001-3103] Regular Pack**

2 x 3 ml syr. and 10 disposable tips

**DIA-ROOT™**
BIO Sealer

BioCermaic Root Canal Sealer

BUY 1, GET 1 FREE!

- Excellent biocompatibility
- Highly antibacterial with high pH 12
- Superior bonding strength to dentin & GP
- Zero shrinkage and low solubility (0.95%)
- Washout resistant to prevent micro-leakage
- Optimal flowability of 22mm
- Setting time: 6 hours
- Easy to store with no wasted material

[1003-701] Dia-Root Bio Sealer Regular Kit
1 x 2g syringe, 20 tips, 1 cap**DiaPaste™**

Calcium Hydroxide with Barium Sulfate

BUY 1, GET 1 FREE!

- 1 Water soluble: easy to remove
- 2 Iodoform-free and biocompatible
- 3 Pre-mixed and pre-filled: save time!
- 4 Antibacterial & bacteriostatic properties
- 5 Highly radiopaque
- 6 Stable: does not solidify nor separate

[A1001-401] Regular Kit 1 x 2g syringe and 20 tips**DiaPex® Plus**

Calcium Hydroxide with Iodoform

BUY 1, GET 1 FREE!

- 1 Pre-mixed paste eliminates mixing
- 2 Pharmaceutical grade iodoform
- 3 Excellent accessibility to apex
- 4 Antibacterial & bacteriostatic properties
- 5 Highly radiopaque
- 6 Does not harden, so $\text{Ca}(\text{OH})_2$ remains active

[A1001-501] Regular Kit 1 x 2g syringe and 20 tips**DUO-GUN™**

Cordless Backfill Obturation Device

**BUY & GET:**
1 BOTTLE GP OBTURATORS
1 DIA-ROOT BIO SEALER**[4007-2001] Duo-Gun Regular Kit****DIA-DUO™**

Complete Cordless Obturation System

**BUY & GET:**
1 BOTTLE GP OBTURATORS
1 BOX DUO-GUN TIPS (23/25G)
1 DIA-ROOT BIO SEALER**[4007-2000] Dia-Duo Regular Kit****DUO-PEN™**

Cordless Warm Vertical Compaction Device

**BUY & GET:**
1 DUO-PEN CONDENSER
1 DIA-ROOT BIO SEALER**[4007-2002] Duo-Pen Regular Kit****Duo-Gun™ Features**

- Variable temperature settings: 140°C, 170°C, 200°C
- Temperature reaches 200°C within 15 seconds
- Designed for both left-handed and right-handed individuals
- Operate hand trigger with minimal effort
- Strong and durable tips are available in both 23G (#4007-2120) and 25G (#4007-2121)

**Duo-Pen™ Features**

- Tightly compacts and seals all canals incl. lateral canals
- 3 temperature settings: 170°C, 200°C, and 230°C
- Softens, spreads, cuts, and compacts gutta percha
- Quick heating tip reaches 230°C within 1 second
- 5 color-coded tip sizes: #.04-40, .04-45, .04-50, .06-50, .08-55