

## D.LUX+

Dual-Wavelength Cordless Curing Light

**BUY 1, GET  
20 DIAFIL OR 40 DIAFIL FLOW FREE\***

Scan for video



- 1 Powerful 2,400 mW/cm<sup>2</sup> output
- 2 Dual wavelength: 385 - 515 nm
- 3 360° swivel light probe
- 4 5 versatile curing modes
- 5 Built-in Light Intensity Sensor
- 6 High capacity 2,600 mAh battery
- 7 Wireless inductive charger

### [4008-1001] D-Lux+ Regular Kit

Contains: D-Lux+ hp & charger, light probe and protector, disposable sheaths (200/box), power supply (adapter)



## DIA-ROOT™ BIO Sealer

**BUY 1, GET 1 FREE!**



- Excellent biocompatibility
- Highly antibacterial with high pH 12
- Superior bonding strength to dentin & GP
- Zero shrinkage and low solubility (0.95%)
- Washout resistant to prevent micro-leakage
- Optimal flowability of 22mm
- Setting time: 6 hours
- Easy to store with no wasted material

### [1003-701] Dia-Root Bio Sealer Regular Kit

1 x 2g syringe, 20 tips, 1 cap



## DIAFIL™

Light-Cured Nanohybrid Composite

**BUY 2, GET 1 FREE OR  
BUY 3, GET 3 FREE!**

\*of the same value

- 1 High tensile strength & low shrinkage
- 2 Easy to place/sculpt and does not slump
- 3 Does not stick to instruments
- 4 Adapts well to margins
- 5 Wear & high fracture resistant
- 6 12 shades with remarkable color stability
- 7 Durable with long-lasting gloss shine



### SYRINGE REFILL (1 x 4g)

Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20, A30



### CAPSULE REFILL (20 x 0.25g)

Shades: A1, A2, A3, A3.5, B2

## DIAFIL™ FLOW

Light-Cured Flowable Composite

**BUY 2, GET 1 FREE OR  
BUY 3, GET 3 FREE!**

- Excellent color adaptability
- Pit & fissure sealant
- Excellent tensile strength and low shrinkage
- Highly radiopaque

### REFILL PACK (2g syringe & 10 tips)

Shades: A1, A2, A3, A3.5, B2, B3, C2, A20, A30.



## DIA-X BOND™

Universal

One-Bottle Universal Bonding Agent

**BUY 1, GET 1 FREE!**

- 1 Versatile with all etching techniques
- 2 All direct and indirect restorations
- 3 Light-, self- & dual-cure compatible
- 4 Virtually no post-op sensitivity
- 5 Advanced 10-MDP monomer
- 6 Convenient: true 1-bottle bond
- 7 Uncompromising bond strength

### [A2001-2201] Dia-X Bond Universal 1 x 5 ml bottle

