




# SciCan OPTIM 1

ONE-STEP CLEANER & DISINFECTANT

Active Ingredient: **0.5% Hydrogen Peroxide**



**1 MINUTE**  
CONTACT TIME

## Fast and Gentle

With **OPTIM 1** there is no need to compromise between efficacy and being protected. With a 1-minute contact time it is one of the fastest and most effective cleaner & disinfectants. It is also safe and gentle for you and surfaces. **OPTIM 1** does not require any GHS hazard pictograms, signal words or statements. **OPTIM 1** provides you with confidence that you are compliant and protecting your patients, staff, equipment and the environment; the perfect balance.



### Effective One-Step Cleaner

- › Removes organic debris while disinfecting, saving you time and money<sup>1</sup>



### Fast Broad-Spectrum Disinfectant

- › 1 minute tuberculocidal, bactericidal, virucidal, fungicidal contact time
- › 1-minute contact time against COVID-19
- › Demonstrated 1-minute virucidal efficacy extends to new and emerging pathogens
- › Ensures compliance, saving you time and money<sup>2</sup>



### Safe for yourself and your patients

- › Non-irritating to skin, eyes and respiratory tract
- › Non-toxic to reproduction, child development and fertility
- › No VOCs (volatile organic compounds)



### Environmentally Friendly

- › Active ingredient, Hydrogen Peroxide, breaks down into water and oxygen
- › No active residues, reducing environmental impact



### Time and Cost Efficient

- › ½ as many wipes as the leading competitor are needed to keep your surface wet long enough to achieve disinfection<sup>3</sup>



### Compatible

- › Designed to be compatible with a wide range of materials<sup>4</sup>

<sup>1</sup> John A. Molinari, PhD, Peri Nelson, BS, and Devi JV. Raghavan, BDS. OPTIM 1 Environmental Surface Cleaning Investigation. Ann Arbor: Dental Advisor, 2016.  
<sup>2</sup> John A. Molinari, PhD, Peri Nelson, BS, Antonia Molinari. Environmental Surface Wetness Test: Comparison of Disinfectant Wipes. Ann Arbor: Dental Advisor, 2015.  
<sup>3</sup> John A. Molinari, PhD, Peri Nelson, BS, Antonia Molinari. Environmental Surface Wetness Test: Comparison of Disinfectant Wipes. Ann Arbor: Dental Advisor, 2015.  
<sup>4</sup> John A. Molinari, PhD, Peri Nelson, BS. Effect of OPTIM 1 on Tensile Properties of Dental Chair Fabrics. Ann Arbor: Dental Advisor, 2016.

**COLTENE**



## Fast Contact Times

A contact time is the amount of time a surface must stay wet for disinfection to occur. OPTIM 1's contact times are faster which ensures compliance and provides confidence that disinfection has been achieved.

OPTIM 1 has been proven effective against difficult to kill pathogens such as Adeno, Polio, and Norovirus, giving you confidence that your facility is protected against pathogens of concern.

Surface Disinfectant Contact Time*																
Product Name	Active Ingredient	Mycobacteria	Non-enveloped Viruses				Enveloped Viruses					Bacteria				Fungi
		TB	Adeno	Polio	Noro	Rota	RSV	HIV	HBV (Hep.B)	HCV (Hep.C)	Influenza A	MRSA	VRE	E.coli	PSA	
OPTIM 1	Hydrogen Peroxide	1	1	1	30 sec	30 sec	30 sec	30 sec	30 sec	30 sec	30 sec	1	1	1	1	1
OPTIM 33TB	Hydrogen Peroxide	1		1	1	1		1	1	1	1	1	1	1	1	3
Birex SE	Phenolics	10	10					10			10	10		10	10	10
Patterson pdCare	Quats-alcohol	3						2	2	2	2	2	2		3	3
Micro-Kleen 3	Quats-alcohol	3		3	3	3		3	3	3	3	3	3		3	3
CaviWipes	Quats-alcohol	3						3		3	3	3	3	3	3	3
Sani-Cloth Super	Quats-alcohol	2	2			2	2	2	2	2	2	2	2	2	2	2
Germicide 3	Quats-alcohol	3		3				3	3	3	3	3	3		3	3
BioText	Alcohol Anhydrous	5													3	3
BioSurf	Alcohol Anhydrous	1													20 sec	10 sec
tb Minuteman	Alcohol Anhydrous	1		1				1			1	1	1		1	1
CaviWipes HP	Hydrogen Peroxide	1	1		1	1		1	1	1	1	1	1	1	1	1
Clorox HP Wipes	Hydrogen Peroxide	1			1	1	30 sec	30 sec	30 sec	30 sec	30 sec	45 sec	30 sec	30 sec	30 sec	4

\*Contact time in minutes. If blank, no published information was available. Information was obtained from publicly available material in January 2024. Information will be updated if there are changes in claims.

Ordering Information for OPTIM 1	
Ready-to-use Pre-saturated Wipes	
SCI1W12	OPTIM 1 RTU Disinfectant Wipes Unscented (12 x 160 sheets)
Liquid	
SCI14X4	OPTIM 1 RTU Disinfectant, Unscented, 4L Bottle (4x4L)
Dry Wipes Refill & Accessories	
SCI33R18	Dry Wipes Refill, 6" x 7" (18 x 160 sheets / case)
SCI33R18XL	Dry Wipes Refill, Large, 10" x 10" (18 x 60 sheets / case)

**For more information contact:**

**SciCan Ltd.**

1440 Don Mills Road  
Toronto, Ontario  
M3B 3P9 / Canada  
T +1 416-445-1600  
F +1 416-445-2727  
TF +1 800-667-7733  
customerservice@scican.com

[www.scican.com](http://www.scican.com)

OPTIM is a registered trademark of SciCan Ltd.

Manufactured by: Virox Technologies Inc. | 2770 Coventry Rd, Oakville, ON L6H 6R1 | Phone: (905) 813-0110

**COLTENE**