



The all new

OVO³

Created by clinicians,
approved by dentists.

Created by clinicians, approved by dentists.

The OVO3 scanner is built on a foundation of over **1,000,000 scans and real-world clinical experience.**

We leverage insights from over 40 reference dentists and collaborate with them daily. In addition to contributions from our internal team of eight dentists, we've gone beyond industry standards, pushing digital innovation **to deliver tangible and relevant clinical benefits.**

We have reassessed and improved every aspect of scanner usage, making it the only truly clinical scanner on the market.

Engineered from the ground up for unmatched accuracy, cutting-edge technology, and effortless usability, this is more than just a tool—**it's the future of clinical precision and efficiency.**



KEY INFORMATION

- French-American company
- Over 1 million scans
- 100% manufactured in the U.S.
- Over 40 reference dentists worldwide



Scan the code to learn more about 3DISC & what drives us.

Go to: www.3disc.com/our-company/



Perfectly connected

USB-C to USB-C

Connect directly and seamlessly to your PC with the new OVO3. Equipped with a proprietary 2-meter USB-C to USB-C cable, it ensures fast, stable data transfer and a consistently reliable connection for uninterrupted performance.



Customizable Button UI

Improved Hygiene

Start scanning with a single click and operate it with one hand, no keyboard or mouse needed. Customizable button functions let you tailor the process to your workflow.



2024 Sensor Technology

4x Higher Resolution | 10-Bit Architecture

Engineered from the ground up for unmatched accuracy our latest optical system offers 20mm scan depth & 4X higher resolution, doubling sharpness and detail in every direction, resulting in clearer images with finer details and smoother color transitions making diagnostics and treatment planning easier and more accurate.



Incredibly Ergonomic

Less than 200g (including tip)

Building on the class leading ergonomic design of OVO, OVO3 is designed for the palm of your hand and provides unmatched comfort and usability. Developed for daily use, it's the perfect balance of form & functionality.



3DiscClinic | A 100% clinical experience

Vibrant & Natural Colors


With a comprehensive set of digital tools, 3DiscClinic guides the dentist through an optimal impression-taking process, ensuring patient comfort and precision restorations, from simple cases to the most complex. Customizable button

functions let you tailor the process to your workflow, while the redesigned user interface provides a 100% clinical experience. Combined with OVO3's 10-Bit architecture, 3DISC's advanced optical imaging technology delivers vibrant & natural colors.

So Fast. So Fluid. So Fly AI.

Automatic Soft Tissue Detection

Fly AI is a groundbreaking upgrade to 3DiscClinic transforming how data seamlessly comes together. Fly AI is designed to automatically detect soft tissue and foreign objects like gloves.



Quality Map



Measurement Tool



Model Builder



Occlusal Clearance



Cross Section Tool



Undercut Tool

3DiscClinic

Intergrations

3DiscClinic enables seamless communication between clinics and labs with one click, offering flexible export options (STL, OBJ, PLY) and

integration with open CAD/CAM systems like DWOS and exocad. Files can be exported locally or via the 3DiscCloud platform.

exocad

3shape


formlabs

dental wings

Clinux

ALPHABYTE

...and more.




Scan&Tell | A patient that sees, agrees.

What is Scan&Tell?

A revolutionary 3D & 2D interactive experience.


Scan&Tell allows the dentist to present the buccal health status & engage with the patient, supporting the dialog in treatment planning.

The Scan&Tell platform unlocks the potential of 1-2 additional treatments per patient.*
Transforming dental care with Scan&Tell – where seamless communication results in increased trust & treatment acceptance.



	Not Using Scan&Tell	Using Scan&Tell
Average number of clinical issues detected	2,4	2,8
Number of accepted treatments	1,3	2,5

- Average number of clinical issues detected
- Number of accepted treatments



	Not Using Scan&Tell	Using Scan&Tell
Crowns & Bridges	6	26
Fillings	23	26
Implants	7	15
Other types of restorations	31	60

- Average number of clinical issues detected
- Number of accepted treatments

1 A statistically significant result was determined with a p-value less than 0.05.
2 Scan&Tell significantly increases the number of crowns, bridges, and implants by +330% (p=0.05) and +114% (p=0.05), respectively.



Scan the QR code & download the download the complete study.

3disc.com/scan-and-tell

“We strive to *give the best user experience*”

Group Customer Experience

3DISC Support

Our customer support team, available by phone or email can access your acquisition PC remotely for a more efficient and effective support experience.

No Annual Fees

A one-time, all-inclusive payment. No additional licensing or case fees, plus FREE software updates.

3DiscCloud

- Streamlined cases submission
- Sleek overview of cases and follow up
- Interactive messaging and discussion
- Available for all users and partners (labs)

Proactive Configuration

OVO3 is custom configured to provide an enhanced user experience and continuous optimization.

3DISC Support

E-Mail: support@3disc.com

From 9 am to 6:00 pm CET
Phone: +33 1 86 65 68 93

From 9 am to 6:00 pm UTC-4
Phone: +1 800 570 0363



Learning Academy

Webinars, Clinical Cases & Tutorials

Develop your knowledge about 3DISC products, the software and all its tools and possibilities. Discover the 3DISC training videos & Webinars.




Learn more about OVO3.

Scan the QR code to access the 3DISC Learning Academy.

3DISC, a global company.

3DISC is a privately-owned company with offices in the United States and France. Our dedicated team consists of digital experts from a wide range of high-tech fields, all working together towards a single goal:

To revolutionize the future of digital dentistry. Our primary goal is to go beyond the conventional limits of digital technology by providing inclusive solutions that offer significant clinical advantages, with a particular focus on intraoral scanners.



3DISC Taking Digital Beyond



Dental
Solutions

1287 Matheson Blvd. E,
Mississauga ON, L4W 1R1
www.dfdentalsolutions.ca
1 (866) 740-8829

