


SPECTRAL VISION. INFORMED DECISIONS.

INTRODUCING THE **VELscope Mantis**


VELscope Mantis is a state-of-the-art oral mucosal examination tool to assist the dental professional in performing comprehensive oral examinations. The Mantis integrates both traditional white light visualization, polarized white light visualization and blue-light-induced fluorescence visualization into the same device.

Innovative Multi-Mode System for Versatile Oral Screening + Discovery



Non-Polarized Light

Non-polarized white light visualization is the traditional method for examining tissue, providing a combination of surface reflections that highlight texture and topography, along with color variations that indicate the state of the underlying tissue.



Polarized Light

Polarized white light visualization excels at showing subtle colour tissue changes by improving colour response and eliminating surface reflections that can be distracting and obscure underlying tissue.



Fluorescence Light

The new VELscope Mantis builds upon its predecessors by increasing field of view and enhancing filtering for a brighter, improved image with increased yellow/orange spectral content.



COMING SOON



JOIN THE WAITLIST:



VELSCOPE

DF Dental
Solutions

1 (866) 740-8829
www.dfdentalsolutions.ca