

### Dia-PROT Gutta Percha

Compliment: Dia-X™ Files, ProTaper® Files, EdgeTaper™

**BUY ANY 4 DIA-PROT GP, GET 1 PACK OF DIA-X FILES 25MM ASSORTED! (#571-291)**



#### { COMPATIBILITY CHART }

F1	ML 150-S601
F2	ML 150-S602
F3	ML 150-S603
F4	ML 150-S604
F5	ML 150-S605
F1-F3	ML 150-S691
F4-F5	ML 150-S692
F1-F5	ML 150-S693

(60/box)

<b>Dia-X Files</b> <b>DX</b>	: Access File (19mm)
<b>Dia-X Files</b> <b>D1</b>	: ProTaper® S1
<b>Dia-X Files</b> <b>D2</b>	: ProTaper® S2
<b>Dia-X Files</b> <b>D3</b>	: ProTaper® F1: Dia-ProT <b>F1</b>
<b>Dia-X Files</b> <b>D4</b>	: ProTaper® F2: Dia-ProT <b>F2</b>
<b>Dia-X Files</b> <b>D5</b>	: ProTaper® F3: Dia-ProT <b>F3</b>

## D.LUX+

Dual-Wavelength Cordless Curing Light

**BUY 1, GET 20 DIAFIL\*!**

\*A2 SYRINGE REFILLS (A2001-1202)

- 1 Powerful 2,400 mW/cm<sup>2</sup> output
- 2 Dual wavelength: 385 - 515 nm
- 3 360° swivel light probe
- 4 5 versatile curing modes
- 5 Built-in Light Intensity Sensor
- 6 High capacity 2,600 mAH battery
- 7 Wireless inductive charger

#### [4008-1001] D-Lux+ Regular Kit

Contains: D-Lux+ hp & charger, light probe and protector, disposable sheaths (200/box), power supply (adapter)



## DIA-ROOT™ BIO Sealer

Bioceramic Root Canal Sealer

**BUY 1, GET 1 FREE!**



**Biocompatible • Antibacterial • Superior Bonding**

#### [1003-701] Dia-Root Bio Sealer Regular Kit

1 x 2g syringe, 20 tips, 1 cap

## DIA-PROSEAL

Resin-Based Root Canal Sealer

**BUY 1, GET 1 FREE\*!**

\*of the same value.



- Excellent sealing ability
- Low shrinkage & solubility
- Fast preparation with dual syringe
- Highly antimicrobial with Ca(OH)<sub>2</sub>

#### [1003-201] Regular Kit

16 g dual-syringe, Mixing Pad and Spatula

#### [1003-202] Intro Kit

4 g dual-syringe, Mixing Pad and Spatula

## e-Temp™

SELF-CURED Temporary Filling Material

- Adheres to dentin
- Expands slightly during setting
- Contains zinc oxide and calcium sulfate
- Excellent adaptability and easy to remove
- Improved mechanical properties for quick temp filling
- Moisture-activated: sets within 25 minutes or less

**BUY 2, GET 1 FREE!**



[#2004-1001] e-Temp (30g) White

## DiaTemp™ DiaTemp™ Flow

LIGHT-CURED Temporary Filling Material

**BUY 2, GET 1 FREE\*!**

\*mix & match okay.

- 1 DiaTemp: high-visc. DiaTemp Flow: low-visc.
- 2 Excellent adhesion
- 3 Great for implant access holes
- 4 Strong, durable, & odorless
- 5 Simple and fast removal in one piece

#### [A2003-1101] BLUE DiaTemp

[A2003-1201] YELLOW DiaTemp

Contains: 3 x 3 g syringe

#### [A2003-2101] BLUE DiaTemp Flow

[A2003-2201] YELLOW DiaTemp Flow

Contains: 5 x 1.8 g syringe and 10 tips



**BUY 2 DIA-X SIL HB OR LB (2 x 50 ml only),  
GET 1 MORE FREE!** (Choice of HB or LB 2 x 50 ml)

Eligible Codes: A2007-1003 & A2008-1003

# DIA-X SIL™

Hydrophilic VPS Impression Material

**NEW!**



### Heavy Body / Light Body

- 1 Excellent elastic recovery
- 2 Hydrophilic properties for excellent flow
- 3 High tear strength from deformation
- 4 Long-lasting dimensional stability
- 5 Excellent detail reproduction
- 6 Patient comfort: easy to remove

[A2007-1003] Heavy Body (2 x 50 ml)  
[A2007-1007] Heavy Body (10 x 50 ml)  
[A2008-1003] Light Body (2 x 50 ml)  
[A2008-1007] Light Body (10 x 50 ml)

### BITE

- 1 Rapid setting time of 1:00 min!
- 2 High hardness: precision & distortion-free
- 3 Patient comfort
- 4 Easy handling
- 5 Outstanding occlusal details
- 6 Clean trimming and grinding

[A2006-1005] Bite Registration  
(4 x 50 ml + 10 mixing tips)



# DIAPLUS™

One-Step 5th Gen Bonding Agent

**BUY 1, GET 1 FREE!**

- 1 Low film thickness
- 2 Direct & indirect applications (5<sup>th</sup> gen)
- 3 No adhesive show-through
- 4 Virtually no post-op sensitivity
- 5 Ethanol-based
- 6 Compatible with dual- & self-cured resin

[A2001-2101] DiaPlus Bond (5ml)

# DIA-X BOND™

One-Bottle Universal Bonding Agent

**Universal**

**BUY 1, GET 1 FREE!**

- 1 Versatile with all etching techniques
- 2 All direct and indirect restorations
- 3 Light-, self- & dual-cure compatible
- 4 Virtually no post-op sensitivity
- 5 Advanced 10-MDP monomer
- 6 Convenient: true 1-bottle bond

[A2001-2201] Dia-X Bond Universal (5ml)



# DIAETCH™

37% Phosphoric Acid Etching Gel

**BUY 2, GET 1 FREE!**

- Fast! 15-second etch
- Precise! 25 gauge pre-bent tips
- Washes off easily with no residue
- Beautiful blue color
- Ideal consistency and biocompatible

[A2001-3103] Regular Pack (2 x 3ml syr. & 10 tips)



# DIAFIL FLOW™

Light-Cured Flowable Composite

**BUY 2, GET 1 FREE!  
BUY 3, GET 3 FREE!**

- Excellent color adaptability
- Pit & fissure sealant
- Excellent tensile strength and low shrinkage
- Highly radiopaque

REFILL PACK (2g syringe & 10 tips)

Shades: A1, A2, A3, A3.5, B2, B3, C2, A20, A30.

